



23910 N 19th Avenue
Building 4 Suite 66
Phoenix, Arizona 85085 USA
1 (800) 243-0820

CLEANING OF RADSCAN SLICKER COVERS AND CUSHIONS

The Slicker is not sterile when shipped. We suggest a thorough cleaning with a standard cleaner or disinfectant commonly used in your hospital environment, or a bleach/water solution as recommended by the Center for Disease Control before each use. Germicide may also be used, however it may cause discoloration and create a sticky film on the cover.

Always leave the CT Slicker or Slicker Cushion on the cradle to maximize table cleanliness.

We suggest a one-year replacement cycle under normal use to provide the highest level of cleanliness for your patients. Under heavy use, the product should be replaced more frequently.

CDC Bleach Disinfection Guideline Links

Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008 (pg. 39-42):
<https://www.cdc.gov/infectioncontrol/pdf/guidelines/disinfection-guidelines-H.pdf>

Guidelines for Environmental Infection Control in Health-Care Facilities:
<https://www.cdc.gov/infectioncontrol/pdf/guidelines/environmental-guidelines-P.pdf>

Manufacturer's brochure on bleach use for disinfection:
<http://www.cloroxprofessional.com/cdiff/pdf/TheTruthAboutBleach.pdf>